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
1.- NAME AND INTENDED USE

NatQal HBV is an HBV DNA-positive control sample. It is formulated for use with *in vitro* qualitative assays, to detect the presence of Human Hepatitis B Virus (HBV) DNA in human plasma.

This device is specifically designed for Nucleic Acid Testing (NAT) assays in blood screening laboratories. NatQal HBV is intended as a run control for routine NAT assays. It allows the user to detect procedural errors or possible failures of the analytical process.

NatQal HBV must be used in accordance with the laboratory routine procedures and reagents manufacturers' recommendations. It must not be substituted for the mandatory positive and negative control reagents provided with manufactured test kits.

NatQal HBV is currently used for NAT assays in blood screening laboratories and participates in the validation of the analytical process. NatQal HBV should be included in every series of NAT assays for HBV DNA. A series of assays (or run) is defined as the number of tests set up at the same time, under the same conditions and processed in a similar manner.

NatQal HBV has been validated on analytical systems supplied by Grifols (Procleix® Tigris® System - Procleix® Ultrio® Assay, Procleix® Panther® System - Procleix® Ultrio Elite® Assay and Procleix® UltrioPlexE Assay), as well as on systems supplied by Roche Diagnostics (cobas® s201 system - cobas® TaqScreen MPX Test, cobas® 6800/8800 systems – cobas® MPX)¹.

2.- SUMMARY

Routine use of independent positive controls enables laboratories to monitor day-to-day and operator variation, lot-to-lot performance of test kits, and can assist the user in identifying errors and failures, and to reveal possible defects in the analytical chain or in the reagents used.

On NAT platforms of French blood screening laboratories (Etablissement Français du Sang (EFS)), NatQal HBV simulates a positive sample or a positive sample pool. It follows all the steps of the analytical process going from the preparation of the sample, to its amplification and detection.

The viral load of the device, evaluated (as an indication) at 30 IU/ml, represents approximately threefold over the 95% detection limit (SD95%) for NAT assays, as determined by the EFS blood screening laboratories¹.

In this context, NatQal HBV participates in the validation of a given series of assays, as well as for the validation on receipt of the reagents used.

3.- REAGENTS

Positive control sample infected with HBV DNA, manufactured from:

- ▶ HBV DNA positive human plasma or serum non-reactive for antibodies to HIV 1 and 2, HTLV I and II, and HCV
- ▶ Pooled human plasma sourced from blood donations collected, qualified and prepared in accordance with good transfusion practices², filtered and tested non-reactive for HIV-1 and HCV RNA and HBV DNA.

It can contain stabilizers and sodium azide as preservative at concentration less than 0.09%

4.- WARNINGS AND PRECAUTIONS

CAUTION

Handle NatQal HBV and all human blood derivatives as though capable of transmitting infectious agents according to Good Laboratory Practice.

For In Vitro Diagnostic Use

Only qualified personnel to perform qualitative and/or quantitative NAT assays should use this device. It must be handled in appropriate containment facilities by competent personnel who are fully trained in the handling of infectious products and in accordance with national and local safety rules.

4.1.- Safety precautions

Rules for the handling of contaminated human blood products must be applied³. Use the World Health Organization (WHO) recommended universal precautions for handling NatQal HBV and human blood.

Do not smoke, eat or drink in areas where specimens are being handled. Do not pipette by mouth. Wear appropriate laboratory protections (lab coats, gloves, goggles, etc...).

Clean any spillage by immediately wiping up with 0.5% sodium hypochlorite solution. After use, dispose of all specimens, controls and materials used in testing as though they contain infectious agents, in accordance with laboratory procedures and Good Laboratory Practice.

4.2.- Handling precautions

Do not use NatQal HBV beyond the expiration date.

Avoid contaminations when opening and closing vials.

To prevent formation of potentially explosive compounds due to reaction of sodium azide and copper or lead pipes, flush waste lines with large quantities of water

5.- PACKAGING

Reference **190015** : NatQal HBV - 1 ml x 50 vials packaging

Reference **190016** : NatQal HBV - 1 ml x 10 vials packaging

The product barcode, which allows computer based reading of the batch code, is included inside the packaging.

6.- PROCEDURE

6.1.- Material required but not provided

Refer to user manuals provided by the manufacturers of NAT kits used.

6.2.- Instructions for use

NatQal HBV must be **quickly and completely thawed at 37 °C** to avoid the formation of cryoprecipitates. Gently mix the tube at this stage until ice clot has disappeared.

Briefly vortex the tube before use. If necessary, remove residual drops of solution in the stopper and open the vial carefully. If necessary, remove bubbles or foam using a sterile pipette. **Do NOT centrifuge the tube.**

NatQal HBV is for single use, and should NOT be reused.

After thawing, the tube should be used immediately.

NatQal HBV should be included in any test run in accordance with the procedure provided by the manufacturer for unknown samples. It must **NOT** be substituted for the control reagents provided with manufactured test kits.

Within the context of the use of these products on the analytical chains of technical platforms of French blood screening laboratories, refer also to the current recommendations defined by EFS in agreement with the reagents manufacturers.

This device is ready-to-use. It has been designed to be integrated directly into the analytical chain at the beginning of the sample preparation step. It does NOT undergo the pooling step. **DO NOT** transfer the product to another recipient.

NatQal HBV should be used in every given series of analyses.

7.- STORAGE INSTRUCTIONS

Store NatQal HBV at -20 °C or colder until use or expiration date.

To prevent leakage, store vials upright.

Indications of reagent instability or deterioration

NatQal HBV is supplied in solid form (storage and transport at a temperature of -20°C or colder). After thawing and homogenisation, it appears as a translucent yellow liquid.

Alterations in physical appearance may indicate instability or deterioration of NatQal HBV.

Solutions that are visibly turbid should be discarded.

8.- NATQAL HBV LIMITATIONS

NatQal HBV must be stored and used in accordance with the instructions given in this document. Test procedures and interpretation of results provided by manufacturers of test kits must be followed closely. Deviations from procedures recommended by test kit manufacturers may produce unreliable results.

The performance characteristics of NatQal HBV have been established only for the detection of HBV

Adverse shipping and storage conditions or the use of outdated products may produce erroneous results.

9.- EXPECTED RESULTS

NatQal HBV should give a **POSITIVE result**, in accordance with experimental data and the manufacturers' instructions for use of NAT reagents for HBV DNA screening, for the interpretation of results for a reactive sample.

If the results are not consistent with those expected, this may indicate an unsatisfactory performance of the analytical process. The user is then advised to proceed to the analysis of the possible causes (operator error, deterioration of test kit reagents, contamination of reagents, faulty performance of equipment, etc...) and to invalidate the entire run if necessary.

10.- SPECIFIC PERFORMANCE CHARACTERISTICS

The viral load of the device is estimated at 30 IU/ml for HBV DNA. As this value is close to the quantification threshold of available kits, and given the variability observed among different manufacturer's assays, different procedures, different lot numbers and different laboratories, the viral load of NatQal HBV is only indicative¹.

The viral load corresponds approximately at thrice the SD95% for NAT techniques.

Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory.

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