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## 1.- NAME AND INTENDED USE

The Medium B19 Monoparametric Run Control for Nucleic Acid Testing (RC DGV MONO B19-M) is a Human Parvovirus B19 (B19V) DNA-positive control sample. It is formulated for use with *in vitro* quantitative assays, to detect the presence of B19V DNA in human plasma.

This device is specifically designed for Nucleic Acid Testing (NAT) assays, and is intended to be used as a run control for routine NAT analytical process. It allows the user to detect procedural errors or possible failures of the analytical process.

This device must be used in accordance with the laboratory routine procedures and reagents manufacturers' recommendations. It must not be substituted for the mandatory control reagents provided with the manufactured test kits.

RC DGV MONO B19-M is currently used for NAT assays in French blood screening laboratories and participates in the validation of the analytical process. RC DGV MONO B19-M should be included in every series of NAT assays for B19V DNA. A series of assays (or run) is defined as the number of tests set up at the same time, under the same conditions and processed in a similar manner.

RC DGV MONO B19-M has been validated on analytical systems supplied by Grifols (PROCLEIX® TIGRIS® and PANTHER® Systems - PROCLEIX® Parvo/HAV Assay), as well as on systems supplied by Roche Diagnostics (cobas® 6800/8800 systems – cobas® DPX).

## 2.- SUMMARY

Routine use of independent positive controls enables laboratories to monitor day-to-day and operator variation, lot-to-lot performance of test kits, and can assist the user in identifying errors and failures, and in revealing possible defects in the analytical chain or in the reagents used.

On NAT platforms of French blood screening laboratories (Etablissement Français du Sang or EFS), RC DGV MONO B19-M simulates a positive sample or a positive sample pool. It follows all the steps of the analytical process going from the preparation of the sample, to its amplification and detection.

The viral load of the device is targeted at 10,000 IU/mL by comparison with the 3<sup>rd</sup> WHO International Standard for B19V DNA<sup>3-4</sup>.

In this context, RC DGV MONO B19-M participates in the validation of a given series of assays, as well as for the validation upon receipt of the reagents used.

## 3.- REAGENTS

Positive control sample infected with B19V DNA, manufactured from :

- ▶ Human plasma positive for B19V DNA calibrated against the 3<sup>rd</sup> WHO International Standard for B19V DNA<sup>3-4</sup>. This plasma is non-reactive for HIV-1, HCV and HAV RNA, HBV DNA, HBs Antigen, anti-HIV-1/2, anti-HCV and anti-HTLV-1/2 Antibodies.
- ▶ Pooled human plasma sourced from blood donations collected, qualified and prepared in accordance with the good transfusion practices<sup>1</sup>, filtered and tested non-reactive for HIV-1, HCV and HAV RNA and HBV and B19V DNA.

## 4.- WARNINGS AND PRECAUTIONS

### For In Vitro Diagnostic Use

**CAUTION**

Handle RC DGV MONO B19-M and all human blood derivatives as though capable of transmitting infectious agents according to Good Laboratory Practices.

Only qualified personnel to perform qualitative and/or quantitative NAT assays should use this device. It must be handled in appropriate containment facilities by competent personnel who are fully trained in the handling of infectious products and in accordance with national and local safety rules.

### 4.1.- Safety precautions

Rules for the handling of contaminated human blood products must be applied<sup>2</sup>. Use the World Health Organization (WHO) recommended universal precautions for handling RC DGV MONO B19-M and human blood.

Do not smoke, eat or drink in areas where specimens are being handled. Do not pipette by mouth. Wear appropriate laboratory protections (lab coats, gloves, goggles, etc...).

Clean any spillage by immediately wiping up with 0.5% sodium hypochlorite solution.

After use, dispose of all specimens, controls and materials used in testing as though they contain infectious agents, in accordance with laboratory procedures and Good Laboratory Practices.

### 4.2.- Handling precautions

Do not use RC DGV MONO B19-M beyond the expiration date.

Avoid contaminations when opening and closing vials

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Reference **B RC 6M-50**: RC DGV MONO B19-M - 1 ml x 50 vials packaging

Reference **B RC 6M-10**: RC DGV MONO B19-M - 1 ml x 10 vials packaging

The product barcode, which allows computer based reading of the batch code, is included inside the packaging.

## 6.- PROCEDURE

### 6.1.- Material required but not provided

Refer to user manuals provided by the manufacturers of the NAT kits used.

### 6.2.- Instructions for use

RC DGV MONO B19-M **must be quickly and completely thawed at 37 °C** to avoid the formation of cryoprecipitates. Gently mix the tube at this stage until ice clot has disappeared.

**Briefly vortex the tube before use.** If necessary, remove residual drops of solution in the stopper and open the vial carefully. If necessary, remove bubbles or foam using a sterile pipette.

**Do NOT centrifuge the tube.**

RC DGV MONO B19-M is for single use, and should NOT be reused.

After thawing, the vial should be used immediately.

RC DGV MONO B19-M should be included in any test run in accordance with the procedure provided by the reagents manufacturers for unknown samples. It must **NOT** be substituted for the control reagents provided with manufactured test kits.

Within the context of the use of these products on the analytical chains of technical platforms of French blood screening laboratories, refer also to the current recommendations defined by EFS in agreement with the reagents manufacturers.

**This device is ready-to-use.** It has been designed to be integrated directly into the analytical chain at the beginning of the sample preparation step or after the pooling step if it is processed.

It does NOT undergo the pooling step. DO NOT transfer the product to another recipient.

RC DGV MONO B19-M should be used in every given series of analyses.

## 7.- STORAGE INSTRUCTIONS

Store RC DGV MONO B19-M at -20 °C or colder until use or expiration date.

To prevent leakage, store vials of RC DGV MONO B19-M upright.

### Indications of reagent instability or deterioration

RC DGV MONO B19-M is supplied in solid form (storage and transport at a temperature of -20°C or colder). After thawing and homogenization, it appears as a translucent yellow liquid.

Alterations in physical appearance may indicate instability or deterioration of RC DGV MONO B19-M. Solutions that are visibly turbid should be discarded

## 8.- RC DGV MONO B19-M LIMITATIONS

RC DGV MONO B19-M must be stored and used in accordance with the instructions given in this document. Test procedures and interpretation of results provided by test kits manufacturers must be followed closely. Deviations from procedures recommended by test kits manufacturers may produce unreliable results.

The performance characteristics of RC DGV MONO B19-M have been established only for the detection of B19V DNA.

Adverse shipping and storage conditions or the use of outdated products may produce erroneous results.

## 9.- EXPECTED RESULTS

RC DGV MONO B19-M should give a **POSITIVE result**, with a viral load targeted at 10,000 IU/mL ( $4 \pm 0.3 \text{ Log}_{10} \text{ IU/mL}$ ) by comparison with the 3<sup>rd</sup> WHO International Standard for B19V DNA<sup>3-4</sup>.

If the results are not consistent with those expected, this may indicate an unsatisfactory performance of the analytical process. The user is then advised to proceed to the analysis of the possible causes (operator error, deterioration of test kit reagents, contamination of reagents, faulty performance of equipment, etc...) and to invalidate the entire run if necessary

## 10.- SPECIFIC PERFORMANCE CHARACTERISTICS

The viral load of the device is targeted at 10,000 IU/mL for B19V DNA, or 4  $\text{Log}_{10} \text{ IU/mL}$ . As an indication, the individual values should be between 3.83 and 4.71  $\text{Log}_{10} \text{ IU/mL}$  ( $\text{CI}_{99\%}$ ).

Procedures for implementing a quality assurance program and monitoring test performance on a routine basis must be established by each individual laboratory.

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