## RedQal

<b>RedQal Group</b>		<b>190045</b>
<b>RedQal Group G</b>		<b>242623</b>
<b>RedQal IAScreen</b>		<b>190046</b>
<b>RedQal IAScreen G</b>		<b>242624</b>
<b>RedQal IAScreen Neg</b>		<b>225860</b>

### Français

#### INTRODUCTION

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### English

#### INTRODUCTION

Internal quality control is mandatory. Any analytical system for blood typing in the ABO Rhesus Kell systems and any test anti-erythrocyte antibodies or irregular antibodies screening must be controlled using a series of samples analysed under the same conditions as those applied to biological samples and for which the results for the relevant analysis are known and guaranteed.

#### PRESENTATION – COMPOSITION

The different devices differ from the type of sample barcodes to allow reading by the claimed automated techniques. For more details on the proposed barcodes, contact the hotline.

##### > RedQal Group and RedQal Group G

This reagent contains 4 x 5.5 mL tubes of whole blood.

These whole blood samples were chosen according to their ABO group and Rhesus Kell phenotype to obtain the following antigenic distribution, which is indicated on the antigram:

##### ABO system :

1 A1 red blood cell  
1 A2 red blood cell  
1 B red blood cell  
1 O red blood cell

##### Rh system :

1 RH red blood cell: -1, -2, -3, 4, 5  
1 RH red blood cell: 1, 2, -3, -4, 5  
1 RH red blood cell: 1, -2, 3, 4, -5  
1 RH red blood cell: 1, 2, 3, 4, 5

At least one of these red blood cells has KEL:1,2 (K+k+) phenotype and another has KEL: - 1,2 (K-k+) phenotype. The KEL:2 (k+) phenotype is deduced for red blood cells with KEL: -1 (K-) phenotype.

##### > RedQal IAScreen et RedQal IAScreen G

This reagent contains 2 tubes of 5.5mL of a plasma or serum dilution with an anti-RH1 antibody and 2 tubes of 5.5mL of a plasma or serum dilution with an anti-FY1 antibody.

Pursuant to the French decree of 15th May 2018, which sets forth the conditions for carrying out erythrocyte immunohematology medical biology examinations, the two samples contain antibodies with a titre less than or equal to 4 in the indirect Coombs filtration technique on a red blood cell having heterozygous expression of corresponding antigen.

For information, the concentration of the anti-RH1 antibody is adjusted to 10 ng/mL\*.

##### > RedQal IAScreen Neg

This reagent consists of two 5.5mL tubes of a plasma or serum dilution free of irregular anti-erythrocyte antibodies.

\* for more details (value for a batch, average dosages, etc.) please contact the UPR Hotline.

#### CONSERVATION – VALIDITY

These samples must be stored with a cap at a temperature between +2°C and +8°C.

Their use must respect the expiry date indicated on each label.

After opening, they must be stored between +2°C and +8°C between each use and until expiry.

Do not freeze.

The presence of a flaky, fibrin-like substance in CQI ABO-D - Rh.K tubes does not interfere with the results. The use of these tubes remains possible.

#### PRECAUTION FOR USE

These reagents are prepared from whole blood collected from voluntary blood donors in accordance with ethical and technical regulations. They were screened for transmissible diseases in accordance with the standards in force set forth in French law. However, despite the negative result of these tests, these samples, as they are biological products of human origin, must still be considered as potentially contaminants and therefore must be handled with the usual precautions.

The use of these devices is strictly reserved for healthcare professionals and more specifically for laboratory technicians who are trained and authorised to perform internal quality control analyses.

Do not use if one or more of the tubes are broken. Waste disposal must be carried out in accordance with current legislation (French decree of 7<sup>th</sup> September 1999 published in the Official Journal of 3rd October 1999).

#### OPERATING INSTRUCTIONS

These samples, identified by a clear number with a bar code, must be used without secondary modifications, under the same analytical conditions as those applied to biological samples in both manual and automatic techniques. Their purpose is to control the analytical process (material - technique - reagent) used for ABO grouping and Rhesus Kell phenotyping and for IAT.

#### EXPECTED RESULTS

For RedQal IAScreen Neg, the result must be negative. For the others devices, the expected results are those clearly written on the antigram or serogram.

#### PERFORMANCE

In recent studies conducted in two EFS user laboratories the RedQal Group was run 288 times and the RedQal IAScreen was run 286 times in screening with untreated panels and 152 times with treated panels. All these tests showed results that were as expected and could be interpreted with no ambiguity.

#### REFERENCES

- French decree of 15<sup>th</sup> May 2018, which set forth the conditions for carrying out erythrocyte immunohematology medical biology tests.

#### SYMBOLS USED

Legal manufacturer

Medical Device for in vitro diagnosis

Batch number

Consult the user manual

Use until

Temperature limits

Catalogue reference

For technical assistance, information or a complaint, please contact us by phone: +33 1-55-93-79-99 or email: [Stde.Clientupr@efs.sante.fr](mailto:Stde.Clientupr@efs.sante.fr)

UPR/PUL/FAB/DF/FI/121  
Version 3– 26/05/2024