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266642

English

INTRODUCTION - APPLICATION OF THE TEST

The Direct Antiglobulin Test is used to demonstrate *in vivo* sensitisation of human red blood cells with anti-IgG and anti-C3d antibody. This test is performed on an anticoagulated sample.

Demonstrating *in vivo* sensitisation relies on the use of human antiglobulin(s) whose Fab portion recognises isotypic markers of immunoglobulins or complement fractions specifically attached to haematopoietic tissue. This test requires the simultaneous and independent use of an IgG antibody and a C3d antibody as well as appropriate reagent controls. The Internal Quality Control for Direct Antiglobulin Testing or IQC DAT is a qualitative test validating the analytical system for the detection of immunoglobulins and/or complement on the surface of red blood cells.

PRESENTATION - COMPOSITION

This reagent consists of two 3.3 mL samples of group O red blood cells suspended at 30% ± 2% in stabilizer medium:

IgG tube: pool of red blood cells sensitised *in vitro* by a human IgG monoclonal antibody

C3d tube: pool of red blood cells sensitised *in vitro* by a complement of human origin

The RedQal Direct and RedQal Direct G differ from the type of sample barcodes to allow reading by the claimed automated techniques. For more details on the proposed barcodes, contact the hotline.

CONSERVATION – VALIDITY

These samples must be stored with a cap at a temperature between +2°C and +8°C, in a vertical position. Their use must respect the expiry date indicated on each label.

After opening, they must be stored at a temperature between +2°C and +8°C until the expiry date.

Slight haemolysis may occur over the time of use without affecting the results. Do not use in cases of marked haemolysis.

Do not freeze nor expose to excessive heat.

PRECAUTION FOR USE

These reagents are prepared from whole blood collected from voluntary blood donors in accordance with ethical and technical regulations. They were screened for transmissible diseases in accordance with the standards in force set forth in French law. However, despite the negative result of these tests, these samples must still be considered as potentially contaminants and therefore must be handled with the usual precautions.

Only personnel authorised to carry out Internal Quality Control analyses should use this device.

Do not use if one or more of the tubes are broken. Waste disposal must be carried out in accordance with current legislation.

OPERATING INSTRUCTIONS

1. RedQal Direct and RedQal Direct G are presented as a 30% ± 2% suspension in stabilizer medium.
2. RedQal Direct and RedQal Direct G are to be used under the same analytical conditions as those applied to biological samples.
3. Upon receipt or prior to first use, tubes must be centrifuged to reform the pack of red blood cells (recommended: 5 min at 2,000 - 3,000 x g). As recommendation: for each new use, the tubes must be homogenised and centrifuged again.

4. For their use, the red blood cells must be suspended in the diluent used for the Direct Antiglobulin Test technique, depending on the support used (refer to the support manufacturer's instructions).

These samples, identified by a clear number with a bar code, must be used without secondary modifications, under the same analytical conditions as those applied to biological samples in both manual and automatic techniques. Their purpose is to control the analytical process (material - technique - reagent) used for the Direct Antiglobulin Test. They were validated using the DC Screening II filtration media from Bio-Rad (Diamed), DG Gel DC Scan from Grifols, Ortho BioVue from Ortho Clinical Diagnostics, and using DC-Lys microplate with the Qwalys technique from Diagast.

EXPECTED RESULTS

		Polyvalent human antiglobulin	Anti-IgG	Anti-C3d	Control
DC Screening II Support	IgG tube	/	Positive	Negative	Negative
	C3d tube	/	Negative	Positive	Negative
DG Gel DC Scan Support	IgG tube	Positive	Positive	Negative	Negative
	C3d tube	Positive	Negative	Positive	Negative
Ortho BioVue anti-IgG / anti-C3b, -C3d / Control Support	IgG tube	/	Positive	≤0.5 +	≤0.5 +
	C3d tube	/	≤0.5 +	Positive	≤0.5 +
DIAGAST DC-Lys Support	IgG tube	/	Positive	Negative	Negative
	C3d tube	/	Negative	Positive	Negative

N.B.: the control must always be conform to the expected results to validate the interpretation.

All other support are excluded from this manual and must be subject to a method validation by the user to define the expected results.

PERFORMANCE

The performance studies were carried out on all the different supports claimed. Inter-assay and inter-batch reproducibility were also assessed. The results demonstrated that the characteristics of the reagents samples allow reproducible performance control in both manual and automated techniques. As part of the batch release process, each batch is tested to guarantee the reagent characteristics described in the table above.

REFERENCES

- French decree of 15th May 2018, which set forth the conditions for carrying out erythrocyte immunohematology medical biology tests.

SYMBOLS USED

Legal manufacturer

Medical Device for *in vitro* diagnosis

Batch number

Consult the user manual

Use until

Temperature limits

Catalogue reference

For technical assistance or request for information, please contact us by phone: +33 1-55-93-79-99 or email: Stde.Clientupr@efs.sante.fr

UPR/PUL/FAB/DC/Fl/001
Version 5 - 26/05/2024