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
---

### INFORMATION FABRICANT / MANUFACTURER INFORMATION

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20, avenue du Stade de France - 93218 LA PLAINE SAINT-DENIS Cedex - Tel. (+33) 6 62 99 46 76 - [efs.reactifs.hotline@efs.sante.fr](mailto:efs.reactifs.hotline@efs.sante.fr)

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For technical assistance or more information, the EFS Reactifs hotline is opened Monday to Friday from 9:00 a.m. to 5:00 p.m.

#### Site de production / Manufacturing site

EFS Hauts de France – Normandie

10/12 Boulevard de Belfort - F-59000 LILLE cedex - Tél : (+33) 3 28 54 22 67 ou (+33) 3 28 54 49 31 - [hfn0.upr-lille@efs.sante.fr](mailto:hfn0.upr-lille@efs.sante.fr)

### 1. INTENDED USE AND INDICATION

SeroQal ALINI 6 is a ready-to-use, unassigned quality control for the markers HBs Ag, HBc Ab, HCV Ab, HIV-1 Ab, HTLV-I Ab and Syphilis Ab. It is available without an assigned value and is ready to use.

It is intended for in vitro use with Abbott ALINITY s reagents and analyzers for serological testing in human plasma / serum. With its low reactivity, it allows for qualitative monitoring of analytical process performance and the detection of procedural errors or potential failures

SeroQal ALINI 6 is an independent control. It must not be substituted for the mandatory control reagents provided with the manufactured test kits.

SeroQal ALINI 6 is an in vitro diagnostic medical device (IVDMD) in accordance with EU Regulation 2017/746.

### 2. COMPOSITION ET PACKAGING

Composition	Description	Packaging	Number of tests
Human plasma and serum collected from blood donors	Positive control for HBs Ag, HBc Ab, HCV Ab, HIV-1 Ab, HTLV-I Ab and Syphilis Ab markers	1 tube 4mL/tube	The sample volume allows at least 1 test for each marker

The human sera and plasma composing SeroQal are controlled with respect to markers HBs Ag, HBc Ab, HCV Ab, HIV p24 Ag, HIV-1 and HIV-2 Ab, HTLV-I and HTLV-II Ab and Syphilis Ab. Samples showing positive serology are treated with methylene blue to inactivate pathogens and thus reduce the risk of infection.

SeroQal are filtered to a 0.2µm (clarifying filtration) and are free of preservatives and antibiotics. They are non sterile.

### 3. DELIVERY AND STORAGE CONDITIONS

SeroQAL must be delivered at a temperature  $\leq -20^{\circ}\text{C}$  (controlled temperature or dry ice transport).

Before opening, SeroQal must be placed and stored upon receipt in an enclosure  $\leq -20^{\circ}\text{C}$ . They are then stable until the expiration date indicated on the label.

After opening, SeroQal can be stored between  $+2^{\circ}\text{C}$  and  $+8^{\circ}\text{C}$  for 7 days or can be refrozen after use without exceeding 3 freeze-thaw cycles. Monitoring the number of cycles is the responsibility of the user.

### 4. PRECAUTIONS AND LIMITS FOR USE

SeroQal ALINI 6 use is strictly reserved for health professionals and more specifically, laboratory technicians trained and authorized to perform analyses on the ALINITY s analysers and to handle of biological samples.

Bad transport, storage or use conditions can lead erroneous results. It is recommended not to use and dispose of products if conditions are not guaranteed.

Do not use if tubes are damaged or if products show physical alteration indicating possible instability or degradation of the product.

### 5. SAFETY INSTRUCTIONS

SeroQal do not contain a hazardous ingredient, no substances with European Community workplace exposure limit values, or substances of very high concern (SVHC) in quantities greater than their respective reporting limit values. Consequently, the product is not supplied with a safety data sheet because the concentrations of chemicals contained in this product do not exceed the thresholds defined in the European Regulation 1907/2006.

Products are composed of blood products of human origin from donors subjected to mandatory blood screening tests according to French regulation. The wearing of personal protective equipment is recommended when handling the product.

SeroQal are not subject to the regulations of Transport of Dangerous Goods with regard to the ADR and IATA except for UN1845 labelling for shipment in dry ice.

SeroQal as well as all media that have been in contact with SeroQal must be handled as potentially infectious products. Disposal conditions and disinfection procedures must comply with current regulations.

### 6. PRINCIPLE OF USE

SeroQal ALINI 6 was designed for use in each series of analyses under the same conditions as those applied to the biological samples

It is formulated to produce a positive reactivity, the signal strength of which may vary depending on the laboratory, analyzer, and reagent lot used. It has no assigned value.

### 7. MATERIALS, REAGENT AND CONSUMABLES

The performance of the SeroQal ALINI 6 is guaranteed with the use of the ALINITY s analytical systems in combination with the materials, reagents and consumables specified in this IFU or validated by the user. The use of a different combination must be validated by the user.

ALINITY s HBsAg - Anti-HBc - Anti-HCV - HIV Ag/Ab Combo - HTLV-I/II - Syphilis

### 8. PROTOCOL

The samples are ready for use. They must be left at room temperature until completely thawed, then **gently homogenized, avoiding foaming**.

The manufacturer's Abbott recommendations for biological samples to be analyzed tests procedure apply to SeroQal ALINI 6.

The tube is labeled with a bar code for computerized lot number reading.

### 9. RESULTS INTERPRETATION

A POSITIVE result is expected: S/Co  $\geq 1,0$  for HBsAg, HBc Ab, HCV Ab, HIV-1 Ab, HTLV-I Ab and SYPHILIS Ab.

If the results are not in accordance with the expected results, it may indicate unsatisfactory performance of the analytical process. It is then advised to implement the non-conformities management procedure of the user laboratory, to analyze the possible causes of this result and to invalidate the series if necessary

### 10. PERFORMANCES

Analytical performance studies were carried out on ALINITYs system in combination with the reagents and consumables identified. Fidelity was determined on 3 batches representative of the batches manufactured: repeatability on 30 values minimum and reproducibility on 30 values minimum. The mean, standard deviation and coefficient of variation were measured. Each batch has been tested on 2 user sites. The results were consistent with the expected results for each marker evaluated and the variability (VC) of less than 10% for each test.

### 11. DETECTION OF ANOMALIES

Any incident related to the device's performance and/or safety must be reported to the manufacturer as soon as possible.